



Effective Biological Cleaning & Hygiene Solutions



Janitorial & Washroom
Surface Cleaning & Hygiene
Hospitality & Catering



Helping you and your business to maintain effective hygiene and infection control standards

Bio Technics are here to help keep you, the people around you and your environment clean and safe from harmful pathogens.

Whether you have a hospitality venue, school, or office, we have a range of high performance products that have been tested extensively.



Contents

Janitorial & Washroom

Bio Clean	Biological multi-purpose surface cleaner	4
Bio TCD	Biological toilet cleaner	5
Bio Blocks	Biological urinal blocks	6
Bio Scale Remover	Biological limescale remover	6
Bio ONA	Biological odour neutraliser	7
Gel Clean	Triple-action floor cleaner	7
Bio Boost	Biological septic tank additive	8
Bio Drain Fresh	Biological drain odour control	8
Bio Treat	Biological portable toilet additive	9
Sobo Gold	General purpose detergent	9
Endurocide Detergent Wipes	General purpose wet wipes	10
Spill Magnet	Super absorbent powder	10

Surface Cleaning & Hygiene

Enduro Hand Sanitiser	Alcohol free hand & skin hygiene	11
Endurocide Virucidal	High level disinfectant	12
AntiBak Tablets	Hard surface disinfectant soluble tablets	13
QFD Cellar Hygiene Spray	Quat free surface cleaner & disinfectant	14
QFD Wipes	Quat free food probe & surface wipes	15
AntiBak Residual	28 day residual surface disinfectant	16

Hospitality & Catering

Grease Control Liquid	Daily biological treatment of Fat Traps	17
Grease Control Pillows	Biological treatment of Fat Traps	18
Grease Control Sachets	Weekly biological treatment of Fat Traps	18

Our Green Credentials

For further information on our range of products, visit bio-technics.com

Application areas

- Suitable for use on all surfaces and all areas, including kitchens and washrooms.
- Toilets & urinals • Sinks & taps
- Showers & baths • Tiles
- Table-tops & work-tops
- Curtains & soft furnishings
- Carpets & linoleum
- Laminate & wood floors
- Ovens & hobs • Mattresses
- Organic stains • Drain odours



Features & benefits

Bio Clean's cutting edge biochemistry offers a new approach to the cleaning of any hard surface.

Containing eco-friendly microorganisms, Bio Clean digests organic debris such as dirt, grime, fats, oils and greases at the source. Odour-causing organics are degraded and stains are removed, leaving the surface deep-clean and with a pleasant, long lasting fragrance.

- Powerful & effective
- Long lasting action
- Environmentally friendly
- Pleasantly fragranced
- Bleach free

How to use

- Dilute at 1 part Bio Clean to 1 part water. For best results, dilute with warm water.
- **To clean hard surfaces:** apply liberally to surfaces, wipe and leave to dry naturally.
- **To remove drain line odours:** apply one squirt nightly to the drain line. For stubborn odours, we recommend the use of Bio Drain Fresh.
- **To remove organic stains on carpets and textiles:** wet surface with warm water, apply and leave to dry.

Note

- The unique blend of bacteria in Bio Clean continues to work long after application. When used regularly, a powerful bacterial action is established that over time will help reduce the amount of Bio Clean you need to use.
- Do not use with bleach, strong acids, caustic / alkaline based products or disinfectants: these will neutralise the bacteria and affect Bio Clean's performance.

Pack sizes

Concentrate	Case size	Order code	Ready to use (50% dilution)	Case size	Order code
5 Litre container	Case of 2	BIO-CLN-5L2	750ml trigger sprays	Case of 12	BIO-CLNR-750TS12
5 Litre container	Case of 4	BIO-CLN-5L4	1 Litre trigger sprays	Case of 12	BIO-CLNR-1TS12
20 Litre container	Sold individually	BIO-CLN-20L	5 Litre container	Case of 2	BIO-CLNR-5L2
200 Litre container	Sold individually	BIO-CLN-200L	5 Litre container	Case of 4	BIO-CLNR-5L4

Application areas

- Designed as a traditional toilet cleaner
- Also suitable for use in urinals



Features & benefits

A bleach free biological approach to the cleaning of toilets.

The natural alternative to bleach, Bio TCD uses friendly bacteria and enzymes to naturally break down any waste present in the toilet bowl.

Bio TCD contains biodegradable cleaning agents and scale removing additives to digest even the most stubborn of stains, leaving your toilet sparkling clean and with a fresh alpine scent.

- Powerful & effective
- Long lasting action
- Environmentally friendly
- Pleasantly fragranced
- Bleach free
- Hazard label free
- Septic tank friendly
- pH 7.5

How to use

Bio TCD not only starts working on contact but continues to do so long after application. Naturally thickened, it remains on the bowl's surface until next flushed

- Apply to the bowl and scrub with toilet brush.
- Leave on surface of bowl for as long as possible to allow biological digestion to take place. Ideally apply last thing in the evening or at the end of the working day.

Note

- The unique blend of bacteria in Bio TCD continues to work long after application. When used regularly, a powerful bacterial action is established that over time will help reduce the amount of Bio TCD you need to use.
- Do not use with bleach, strong acids, caustic / alkaline based products or disinfectants: these will neutralise the bacteria and affect Bio TCDs performance.

Pack sizes

Description	Case size	Order code
750ml applicator bottles	Case of 12	BIO-TCD-750ML12
5 Litre container	Case of 2	BIO-TCD-5L2
5 Litre container	Case of 4	BIO-TCD-5L4
20 Litre container	Sold individually	BIO-TCD-20L

Bio Blocks Biological Urinal Blocks

Biological Cleaning

Application areas

- Designed for use in traditional urinals and urinal troughs
- Also suitable for use in "waterless" systems



Features & benefits

Biological, biodegradable and bleach free de-odourising urinal blocks that naturally break down urine and organic waste whilst bringing a touch of fragrance to your washroom.

- Powerful & effective
- Environmentally friendly
- Bleach free
- Pleasantly fragranced
- Septic tank friendly
- pDCB free

How to use

- Individual urinals: place 1 in each urinal.
- Trough urinals: place approximately 50cm apart.
- Replace when reduced to 5-10% of original size.

Pack size

Description	Block size	Quantity	Order code
Tub of blocks	Ø 38mm x h.15mm	Tub of 50	BIO-BLK-50T

Bio Scale Remover Biological Limescale Remover

Biological Cleaning

Application areas

- Designed to remove limescale from toilet bowls
- Also suitable for use in urinals and drainlines.



Features & benefits

A water-soluble sachet designed to clean toilet bowls and combat limescale build up: Bio Scale Remover's biological formulation 'digests' and removes limescale, organic matter and malodours.

- Powerful & effective
- Environmentally friendly
- Bleach free
- Easy to use
- Septic tank friendly
- Helps keep drain lines clear

How to use

- Add one sachet to the toilet bowl and leave to dissolve
- After around 1 minute, scrub bowl and flush.
- Repeat weekly.

Pack size

Description	Sachet size	Quantity	Order code
Tub of 50g sachets	8 x 8 cm	Tub of 100	BIO-SCR-100

Bio ONA Biological Odour Neutraliser

Biological Cleaning

Application areas

- Effective against odours in washrooms and caused by pets & animals, food and tobacco smoke
- Also suitable for use against flood and fire damage odours



Features & benefits

A biological spray that naturally digests the source of odours from hard surfaces and textiles to permanently remove them, leaving a fresh citrus fragrance.

- Neutralises odours
- Environmentally friendly
- Long lasting biological action
- Hazard label free

How to use

- Dilute at 1 part Bio ONA to 2.5 parts water.
- Use in a trigger spray or via automatic spray systems.
- Apply direct to source and repeat daily until neutralised.

Pack sizes

Description	Pack sizes		
Concentrate	2 x 5 Litre	4 x 5 Litres	20 Litres
Ready to use	12 x 1 Litre trigger sprays	2 x 5 Litre	4 x 5 Litres

For further details on packs sizes and order codes, visit bio-technics.com

Gel Clean Triple Action Floor Cleaner

Biological Cleaning

Application areas

- Suitable for use on floor substrates, including linoleum, vinyl, laminates, washroom floors, tiles and varnished wood.



Features & benefits

Gel Clean is an antimicrobial detergent floor cleaner that leaves your floors not just clean, but sanitised, odour-free and polished.

Part floor cleaner and part sanitiser, Gel Clean combines antimicrobial cleaning and surface polishing agents to leave floors powerfully deep cleaned and gleaming.

- Triple action
- Bleach free
- Fresh pine fragrance

How to use

- Add 60ml to 5 litres of warm water and stir to dissolve.
- Clean floor, leave to dry then lightly buff to restore sheen.

Pack sizes

Pack sizes		
2 x 5 Litre	4 x 5 Litres	20 Litres

For further details on packs sizes and order codes, visit bio-technics.com

Application areas

- Designed as a biological septic tank additive.
- Apply via sink, drain or toilet.



Features & benefits

One of the most effective septic tank additives available!

Bio Boost's powerful formula works in harmony with the tank's natural bacteria to prevent crusting and blockages from forming - helping to minimise odours and to reduce costly system failures. Bio Boost will even biologically degrade waste paper in 4-5 days.

- Powerful & effective
- Prevents crusting & blockages
- Long lasting action

Pack sizes

Description	Case size	Order code
1 Litre Bottle	Case of 12	BIO-BST-1L
5 Litre container	Case of 2	BIO-BST-5L2
5 Litre container	Case of 4	BIO-BST-5L4
20 Litre container	Sold individually	BIO-BST-20L

Note

- For best results, do not use with bleach, strong acids, caustic / alkaline based products or disinfectants when using Bio Boost: these will neutralise the bacteria and affect its biological performance. We recommend the use of our biological surface cleaner [Bio Clean](#) and biological toilet bowl cleaner [Bio TCD](#) to complement Bio Boost.

Application areas

- Designed to remove build up of fats and greases from drain lines, one of the main causes of malodours and blockages in drains.



Features & benefits

An easy to use biological drain odour control solution.

Coating the inside of drain pipes with a highly concentrated bio-film, Bio Drain Fresh captures and digests fats and greases, preventing them from building up and causing odours. Bio Drain Fresh's long lasting biological action helps keep drain lines running free and leaves a refreshing fragrance.

- Powerful & effective
- Biologically digests fats & greases
- Long lasting action
- Environmentally friendly
- Bleach free

Pack sizes

Description	Case size	Order code
1 Litre bottle	Case of 12	BIO-DRF -1L12
5 Litre container	Case of 2	BIO-BST-5L2
5 Litre container	Case of 4	BIO-BST-5L4
20 Litre container	Sold individually	BIO-BST-20L

Application areas

- Designed as a biological portable toilet additive
- Also suitable as a rinse aid



Features & benefits

Bio Treat rapidly degrades urine and waste, replacing odours with a pleasant long lasting fragrance.

An environmentally friendly alternative to formaldehyde-based products, Bio Treat rapidly degrades urine and waste and replaces odours with a pleasant long lasting fragrance.

- Long lasting action
- Reduces emptying cycles
- Formaldehyde free

Pack sizes

Description	Case size	Order code
1 Litre Bottle	Case of 12	BIO-TRT-1L
5 Litre container	Case of 2	BIO-TRT-5L2
5 Litre container	Case of 4	BIO-TRT-5L4
20 Litre container	Sold individually	BIO-TRT-20L

Note

- For best results, do not use with bleach, strong acids, caustic / alkaline based products or disinfectants when using Bio Boost: these will neutralise the bacteria and affect its biological performance. We recommend the use of our biological surface cleaner [Bio Clean](#) and biological toilet bowl cleaner [Bio TCD](#) to complement Bio Boost.

Application areas

- Floors
- Walls
- Equipment
- Hard surfaces



Features & benefits

Sobo Gold is a powerful, biodegradable, multi-purpose cleaner and degreaser, ideal for the removal of oil, grease, dirt, grime, drips, spills and splashes from hard surfaces and equipment.

- Powerful & effective
- Biodegradable
- Low to medium foamer
- VOC free

How to use

- Apply by mop and buckets, cloths and spray, pressure washer, floor scrubbing machinery or other cleaning equipment.
- Allow a contact time of 2-3 minutes, before rinsing with clean water

Pack sizes

Description	Case size	Order code
1 Litre bottle	Case of 12	SOBO-GD-CONC-1L
5 Litre container	Case of 4	SOBO-GD-CONC-5L
5 Litre container	Case of 4	SOBO-GD-CONC-25L
1 Litre READY TO USE container	Case of 12	SOBO-GD-RTU-1L

Endurocide Detergent Wipes

General Cleaning

Application areas

- Door handles
- Wall and floors
- Toilets, sinks and taps
- Ventilation grilles
- Work stations
- Bed frames & mattresses



Features & benefits

Multi-purpose, water-based detergent wipes for the everyday cleaning and damp dusting of surfaces, equipment and furniture.

- Alcohol free
- pH neutral
- Fragrance free
- Chlorine free
- Bleach free
- Hazard label free



Pack sizes

Description	Wipe size	Case Size
Bucket of 225 wipes	240 x 250mm	Sold individually
225 wipes Refill	240 x 250mm	Case of 6

Spill Magnet Super-Absorbent Spill Powder

Spill Control

Application areas

- Ideal for use in catering operations, educational institutions, care homes, healthcare facilities, sports centres, shops and public spaces.



Features & benefits

Spill Magnet quickly absorbs and retains unpleasant and hazardous liquid spills. Ideal for use on spills of blood, vomit and urine, Spill Magnet forms a gel which can be easily swept up for safe disposal.

- Fast acting
- Non-leaching
- Spill containing
- Environmentally friendly

How to use

- Scoop Spill Magnet around and over the spill
- Allow 2-3 minutes for Spill Magnet to expand into a gel
- Brush in gently to produce a dry, sweepable material
- Once removed, clean and disinfect the affected area using a suitable product such as Endurocide Virucidal, see page 12
- Dispose of waste according to risk of material spilt

Pack sizes

Description	Case size
5kg bucket	Individual

Enduro Hand Sanitiser Water-based & alcohol free

Hand Hygiene

Application areas

- Designed for everyday regular use - any time, anywhere!
- Catering & kitchens
- Nurseries & playgroups
- Schools & colleges
- Shopping & eating out
- At work & commuting
- Before & after eating
- Before & after visiting the toilet



Features & benefits

Fast acting, quick drying and long lasting, alcohol free Enduro Hand Sanitiser is designed for everyday use - any time, anywhere!

With 60 second kill rates against bacteria, fungi and enveloped viruses, Enduro Hand Sanitiser is dermatologically tested, suitable for repeated use and has proven antimicrobial activity for up to 4 hours after application.

- Alcohol free
- 60 second kill rates
- Dermatologically tested
- Bactericidal
- Fungicidal
- Virucidal*
- Solvent free
- Fragrance free
- No sticky residue
- Long lasting protection
- Naturally moisturising

Test results

Tested according to International Standards:

Pathogen	Standard	
Enveloped viruses	EN 14476	
Bacteria	EN 14561	EN 1276
Fungi	EN 1650	EN 1275
Spores	EN 13704	

Other results

- Dermatologically tested: 96 hour skin patch test done on 30 volunteers
- Effective in 30 seconds: Tested against bacteria and shown to be highly effective after 30 seconds
- 2 hour residuality: Tested against bacteria, fungi and spores and shown to provide up to 99.99% protection 2 hours after application on surfaces
- 4 hour residuality: Tested against bacteria and shown to provide up to 99.8% protection 4 hours after application on surfaces

Full details on our test data is available to professional users on request.

* Testing against enveloped viruses only.

Pack sizes

Dispenser bottles	Case size	Order code	Dispenser refills	Case size	Order code
50ml foam dispenser	Case of 12	HS-50F	1 Litre bottles	Case of 12	HS-1L_RF12
100ml foam dispenser	Case of 40	HS-100F	5 Litre containers	Case of 2	HS-5L_RF2
150ml foam dispenser	Case of 40	HS-150F	5 Litre containers	Case of 4	HS-5L_RF4
200ml foam dispenser	Case of 32	HS-200F			
500ml foam dispenser	Case of 12	HS-500F	Wall dispensers	Case size	Order code
			Manual dispenser (0.8L pouch fill)	Individual	MANW-800P
			Manual dispenser (1L reservoir fill)	Individual	MANW-1000R

Wall dispenser refills and bulk containers are available: for further details and order codes, visit bio-technics.com

Application areas

- High Level Disinfection applications in non-food areas
- Ideal for sports centres, football stadiums, and public spaces



Features & benefits

Endurocide Virucidal provides quick and easy High Level Disinfection on hard surfaces.

- Fully tested to International Standard EN 14476 against non-enveloped viruses: Adenovirus, Poliovirus & MNV Norovirus
- Fully tested to International Standard EN 14476 against enveloped viruses including Coronavirus
- Tested under dirty conditions against bacteria and fungi

Application Advice

- Available as ready to use liquid or as a concentrate (for dilution at 5% or 1 part to 19 parts water)
- Exhibits excellent cleaning properties - for best results and for High Level Disinfection, pre-cleaning of surfaces is recommended
- Can be applied on surfaces by mop, wipe, lance, spray and surface flooding or by immersion for tools and equipment

PACK SIZES

Ready to use	Case size	Order Code
• 1 Litre trigger spray	12	EV-RTU-1
• 5 Litre refill	2	EV-RTU-2X5
• 5 Litre refill	4	EV-RTU-4X5

Concentrate	Case size	Order Code
• 1 Litre bottle	12	EV-C-1
• 5 Litre container	2	EV-C-2X5
• 5 Litre container	4	EV-C-4X5
• 20 Litre container	Individual	EV-C-200

Tested against:
 ☑ Viruses
 ☑ Bacteria
 ☑ Fungi

Application areas

- General purpose surface disinfection in catering operations, educational institutions, care homes or healthcare facilities
- Front of house catering disinfection
- Disinfection of food trays, counters, cabinets & equipment
- Terminal High Level Disinfection
- Decontamination of small scale hazardous spills
- Legionella decontamination of shower heads



* Screen testing only
 ** Testing done using 2 tablets in 500ml of water

Features & benefits

Chlorine free AntiBak Tablets offer an immediate response disinfection solution, replacing chlorine tablets for quick and easy surface cleaning and sanitising.

Odourless and taint free, each AntiBak Tablet provides 500ml of powerful disinfectant proven to be highly effective against viruses,

- 60 second kill rates
- Easy to use
- Ready in 2 minutes
- Virucidal
- Bactericidal
- Fungicidal
- Chlorine free
- Fragrance free
- Taint free
- Solution remains active for up to 7 days

Test results

Tested according to International Standards:

	Pathogen	Standard	Test Time
Viruses	Adenovirus	EN 14476	5 min.
	Norovirus	EN 14476	1 min.
	Bovine Coronavirus	EN 14476*	30 sec.
	H1N1 Swine Flu	EN 14476*	5 min.
	HIV	Independent	2.5 min.
Bacteria	Hepatitis C**	Independent	10 min.
	<i>Campylobacter jejuni</i>	EN 1276	1 min.
	<i>Enterococcus hirae</i>	EN 1276	1 min.
	<i>Escherichia coli</i>	EN 1276	1 min.
	<i>Klebsiella pneumoniae</i>	EN 1276	1 min.
	<i>Legionella pneumophila</i>	EN 13727	1 min.
	<i>Listeria monocytogenes</i>	EN 1276	1 min.
	Methicillin-resistant <i>S.aureus</i>	EN 1276	1 min.
	<i>Pseudomonas aeruginosa</i>	EN 1276	1 min.
	<i>Salmonella enteritidis</i>	EN 1276	1 min.
Fungi	<i>Aspergillus niger**</i>	EN 1650	15 min.
	<i>Candida albicans</i>	EN 1650	10 min.

Other results

- BS ISO 4120-2004: Tested at Campden & Chorleywood Food Research Association Group using the Triangle Test Method TES-S-004 to determine whether the product has the potential to taint food (pass)
- 7 day active solution: Solution made up and tested daily for seven days using EN 1276 against *E.coli* (pass)

All testing carried out using 1 tablet in 500ml of water unless otherwise stated.

Pack sizes

Description	Case size	Order code
Tub of 100 x 5g tablets	Sold individually	AB-T-100

QFD Cellar Hygiene Spray Cleaner & Disinfectant

Hard Surfaces

Application areas

- Ideal for deep cleaning
- Catering operations, front & back of house
- Kitchens, food trays & probes
- Food processing surfaces
- Sinks and taps
- Keg couplers
- Beer tap nozzles and drip trays
- Bar counters and worktops
- Ice buckets and scoops
- Door handles and light switches
- Shelving and storage racks

Features & benefits

A high-performance hard surface cleaner and disinfectant, ideal for use on all hard surfaces in cellars, kitchens and food preparation areas or anywhere requiring a high level of cleanliness and hygiene.

Odourless, taint-free and Quat free, QFD is highly effective against a broad spectrum of fungi, yeasts, moulds and bacteria.

- One step cleaner & disinfectant
- No rinsing required
- Bactericidal
- Fungicidal
- Virucidal*
- Quat (QAC) free
- Food safe - no Maximum Residue Level
- Taint free
- Fragrance free
- Hazard label free

Test results

Tested according to International Standards:

	Pathogen	Standard	Test Time
Viruses	Bovine Coronavirus	EN 14476	1 min.
Bacteria	<i>Campylobacter jejuni</i>	EN 13697	3 min.
	<i>Enterococcus hirae</i>	EN 1276	5 min.
	<i>Enterococcus hirae</i>	EN 13697	2 min.
	<i>Escherichia coli</i>	EN 1276	30 sec.
	<i>Escherichia coli</i>	EN 13697	2 min.
	<i>Legionella pneumophila</i>	EN 13697	5 min.
	<i>Listeria monocytogenes</i>	EN 13697	2 min.
	<i>Pseudomonas aeruginosa</i>	EN 1276	5 min.
	<i>Pseudomonas aeruginosa</i>	EN 13697	3 min.
	<i>Salmonella choleraesuis</i>	EN 1276	5 min.
	<i>Salmonella enteritidis</i>	EN 13697	2 min.
	<i>Staphylococcus aureus</i>	EN 13697	2 min.
	<i>Staphylococcus aureus</i>	EN 1276	5 min.
	Vancomycin-resistant <i>Enterococcus</i>	EN 13697	2 min.
Fungi	<i>Aspergillus niger</i>	EN 13697	15 min.
	<i>Candida albicans</i>	EN 13697	5 min.

All testing carried out under dirty conditions except EN 14476 done under clean conditions.

Pack sizes

Concentrate	Case size	Order code	Ready to use (10% dilution)	Case size	Order code
5 Litre container	Case of 2	QFDCH-C-5L2	750ml trigger sprays	Case of 12	QFDCH-RTU-750TS12
5 Litre container	Case of 4	QFDCH-C-5L4	1 Litre trigger sprays	Case of 12	QFDCH-RTU-1TS12
20 Litre container	Sold individually	QFDCH-C--20L	5 Litre container	Case of 2	QFDCH-RTU-5L2
200 Litre drum	Sold individually	QFDCH-C-200L	5 Litre container	Case of 4	QFDCH-RTU-5L4
			20 Litre container	Sold individually	QFDCH-RTU-20L
			200 Litre drum	Sold individually	QFDCH-RTU-200L

* Testing against enveloped viruses only.

QFD Wipes Alcohol free & Quat free probe & surface wipes

Hard Surfaces

Application areas

- Designed for cleaning and disinfecting items & areas including:
 - Food probes
 - Cold stores
 - Food trays
 - Tabletops
 - Machinery & equipment



Quat free and alcohol free, QFD Wipes keeps food preparation areas clean and disinfected at all times. Also perfect for use on food probes with no rinsing required.

- Bactericidal
- Fungicidal
- Virucidal*
- No rinsing required
- Quat (QAC) free
- No MRL
- Taint free
- Fragrance free
- Hazard label free

Tested according to International Standards:

Viruses: EN 14476 Bacteria: EN 1276 & 13697 Fungi: EN 13697

For full testing information, please refer to QFD Cellar Hygiene on page 14.

Pack sizes

Description	Wipe size	Case Size	Order code
Tub of 200 wipes	130 x 130mm	12 per case	QFD-W-200T
Bucket of 225 wipes	240 x 250mm	Sold individually	QFD-W-225B



AntiBak Residual 28 Day Residual Surface Disinfectant

Hard Surfaces

Application areas

- High contact & traffic areas
- Washroom furniture, sinks & taps
- Door handles
- Light switches
- Kitchen floors
- Foot mats
- Equipment
- Cold stores
- Difficult to access areas

Note: not suitable for use on food contact surfaces



* Tested against enveloped viruses only
** Screen testing only

Features & benefits

A new approach to surface disinfection: AntiBak Residual not only kills pathogens on contact, but continues to do so for up to 28 days after application.

- Tested and proven to be effective 28 days after application
- Bactericidal
- Fungicidal
- Virucidal*
- >99.99% inactivation of Bovine Coronavirus in 30 seconds
- Suitable for use with spray bottles
- Dilute at 1:4 (20%) for Intermediate level disinfection or 1:19 (5%) for Low level disinfection.

Test results

Tested according to International Standards:

	Pathogen	Dilution	Standard	Test Time
Viruses	Bovine Coronavirus	20%	EN 14476	30 sec.
	H1N1 Swine Flu	25%	EN 14476**	5 min.
Bacteria	<i>Campylobacter jejuni</i>	4%	EN 1276	1 min.
	<i>Enterococcus hirae</i>	4%	EN 1276	5 min.
	<i>Escherichia coli</i>	4%	EN 1276	5 min.
	<i>Klebsiella pneumoniae</i>	4%	EN 1276	1 min.
	<i>Legionella pneumophila</i>	10%	EN 13727	10 min.
	<i>Pseudomonas aeruginosa</i>	20%	EN 1276	1 min.
	<i>Salmonella enteritidis</i>	4%	EN 1276	1 min.
Fungi	<i>Salmonella typhimurium</i>	4%	EN 1276	5 min.
	<i>Staphylococcus aureus</i>	4%	EN 1276	5 min.
	Vancomycin-resistant <i>E.faecalis</i>	20%	EN 1276	1 min.
	<i>Aspergillus niger</i>	20%	EN 1650	15 min.
	<i>Candida albicans</i>	20%	EN 1650	15 min.

Other results

- 28 days residuality: Tested against *E.coli* and *Listeria monocytogenes* and proven to remain effective 28 days after application.
- 7 days residuality: Tested against six types of bacteria and fungi and proven to remain effective after 24 hours, 48 hours and 7 days.

Pack sizes

Concentrate	Case size	Order code	Ready to use (20% dilution)	Case size	Order code
5 Litre container	Case of 2	AB-RES-2X5	1 Litre trigger sprays	Case of 12	AB-RES-RTU1
5 Litre container	Case of 4	AB-RES-4X5	5 Litre container	Case of 2	AB-RES-RTU2X5
20 Litre container	Sold individually	AB-RES-20	5 Litre container	Case of 4	AB-RES-RTU4X5
200 Litre drum	Sold individually	AB-RES-200			
1000 Litre IBC	Sold individually	AB-RES-1000			

Grease Control Liquid Daily Treatment

Biological Grease Control

Application areas

- Commercial Kitchens
- Takeaway Outlets
- Drain Lines and Fat Traps



Features & benefits

Grease Control Liquid digests the difficult to degrade fatty acids and organic wastes found in Grease Traps and Drains.

Perfect for those who want to use green products without a compromise on performance, A safe to use blend of micro-organisms specially selected for their ability to degrade both the short and the long chain fatty acids found in restaurant and kitchen Grease Traps.

- Liquefies fats oils and grease
- Reduces blockages
- Reduces malodours
- Low dosage needed
- Bleach free

How to use

- Apply manually or use with Bio Technics Peristaltic Pump
- For best results, apply 2 hours after the end of a working day

Grease Traps

- Application rates vary with the size of the Grease Trap, dose between 300-1000ml dependent on the capacity of grease trap

Slow Moving Drains (Smaller Operations)

- Pre-treat for 3 days at 70ml per day
- After pre-treat, apply 70ml weekly

Slow Moving Drains (Larger Operations)

- Pre-treat for 3 days at 200ml per day
- After pre-treat, apply 200ml weekly

Pack sizes

Description	Case size	Order code
5 Litre container	Case of 4	GC-LQD-5L
20 Litre container	Sold individually	GC-LQD-20
200 Litre Container	Sold individually	GC-LQD-200

Grease Control Pillows Biological Treatment

Biological Grease Control

Application areas

- Commercial Kitchens
- Takeaway Outlets
- Drain Lines and Fat Traps



Features & benefits

An environmentally safe biological pillow that liquefies grease in fat traps used by smaller catering establishments.

- Removes malodours
- Long lasting biological action
- Environmentally friendly
- Time saving

How to use

- Clean Fat Traps and re-seed using the Grease Control Pillow, replace every 8 to 12 weeks
- Use alongside Grease Control Sashets for increased performance

For further details on packs sizes and order codes, visit bio-technics.com

Grease Control Sachets Biological Treatment

Biological Grease Control

Application areas

- Commercial Kitchens
- Takeaway Outlets
- Drain Lines and Fat Traps



Features & benefits

Designed for ease of use and for the smaller catering application, Grease Control Sachets quickly digests Fats, Oils and Greases (FOG), reducing malodours and so reducing the need to manually clean fat traps..

A blend of powerful bacteria powders selected to digest fat, oils and grease in drain lines and fat traps.

- Biological
- Time saving
- Removes malodours

How to use

- Simply empty the contents of a 200g Grease Control Sachet into a 10 litre bucket of warm water and pour down drain, repeating weekly.

Pack sizes

Description	Units per case	Order code
Tub of 200g sashets	12 per case	GC-SCT-200

For further details on packs sizes and order codes, visit bio-technics.com

Our Green Credentials

Caring for the environment is tightly woven into our DNA.

Not only are our products carefully engineered to be highly effective with minimal impact on the environment, but our operations are also strictly monitored to make sure that we as a company, are as environmentally friendly as we can possibly be.

Our company values are simple:

- 100% product satisfaction - guaranteed
- Market-leading products - which are environmentally friendly
- Products tested and accredited to the highest possible standards
- World class customer service and distribution



Our commitment to effective, environmentally friendly products led to the creation of sporicidal products that are:

- Water-based
- 100% free from chlorine, peracetic, bleach & alcohol
- Tested and approved to international standards



Our premises use electrically powered forklifts to reduce emissions



We recycle 97.2% of our waste

The national average is 46%



Our premises is accredited to Environmental Management System ISO 14001:2015



Our premises is heated using a biomass boiler and lit by energy efficient LED lighting



Manufactured by



Bio Technics Ltd.
Linton Business Park, Gourdon, Aberdeenshire, Scotland, UK DD10 0NH
Tel: +44 (0) 1561 361515 info@bio-technics.com bio-technics.com

Proudly manufactured in the UK for over 40 years



Use biocides safely. Always read the label and product information before use. **Disclaimer:** Every effort has been made to ensure that the information presented is correct at the time of publication. However, Bio Technics Ltd does not assume and hereby disclaims any liability to any part for any loss, damage, or disruption caused by errors or omissions resulting from negligence, accident, or any other cause. Bio Technics Ltd makes no reference to any specific brand of competitor products. The data supplied is for guidance purposes only and is therefore without guarantee. Any opinions expressed in this publication are not necessarily those of Bio Technics Ltd and are not for onward publication without consent.